

Supplementary Figure 2. Nuclear extracts (2 μ g) of SIE cells treated for 20 minutes with 50 ng/mL rDrm were incubated with a biotin-labeled NF-kB double-stranded DNA oligonucleotide probe (5'AGTTGAGGGGACTTTCCCAGGC) in the absence or in the presence of anti-RelA/p65 antibodies or of a molar excess of the unlabeled NF-kB probe or of a mutant NF-kB probe (5'AGTTGAGGCGACTTTCCCAGGC). The protein-DNA complexes were analyzed by EMSA onto a native 6% polyacrylamide gel.